




PA-500A Ventilator

★ Features ★

- ◇ It is a universal, poly-mode and multi-functional unit which is commonly used in clinical first aid and respiratory remedy.
- ◇ Respiratory mode: Controlling Assistant/Controlling, synchronous intermittent mandatory, Manual controlling.
- ◇ It possesses respiratory functions as Positive pressure at the end of ventilation, sigh, etc.
- ◇ It is an essential respiratory unit for hospitals.



★ Specifications ★

- | | |
|--|--|
| ◇ Tidal Volume Adjustment: 50~1400ml | ◇ Air Flow Adjustment: 1~12 times/min |
| ◇ Ventilation Capacity per Minute: ≥18L/m(MV) | ◇ PEEP: 0~1.0kPa |
| ◇ O ₂ Concentration: 40%~90% | ◇ Sigh Volume: 150% of tidal volume |
| ◇ Ventilator Conformance: ≤30mL/kPa | ◇ O ₂ Consumption: ≤1.5mPa change in pressure after 1-hour operation of O ₂ cylinder of 12250kPa/40L |
| ◇ Respiratory Rate: 1~99/min | |
| ◇ I/E: 4.0:1.0, 3.0:1.0, 2.0:1.0, 1.0:1.0, 1.0:1.5, 1.0:2.0, 1.0:2.5, 1:3.0, 1.0:4.0 | |
| ◇ Maximum Safe Pressure: ≤6.0kPa | |
| ◇ Ptr. Range: -0.4~ ± 1.0kPa | |
| ◇ Transfer Time from controlled to Assistant Respiration: 6s | |

PERLONG GROUP



Add: No.28-1 Dajiaochang Road, Nanjing, 210022 China
 Tel: 86-25-52635350 52635360 Fax: 86-25-52635356
 http://www.perlong.com E-mail: nanjingpt@perlong.com.cn