



HEALTH CARE
unicos
ACP-700
Auto Chart Projector



Outstanding feature




• Number 1 Brand in the world
The first developed in Korea by the highest technology and accumulated production experience were qualified by world wide market (Over 40 countries) & global OEM business



• Clear and precise projection chart
Four different types of chart show clear and bright images for your right vision test.



• A comfortable operation
A remote controller and fast chart change make operator comfortable and stable.



Specifications

Chart	31pcs
Chart selection speed	Average 0.2sec.
Mask	1 open
	5 horizontal lines
	5 Vertical lines
	21 single letters
Filter	Red/Green
Projection distance	2.0m~7.0m (5m is standard)
Projection magnification	30 x (at 5m)
Projection Size	φ 252mm and 330W x 225H (at 5m)
Tilt angle	Ball joint (±10°)
Lamp	6V 20W halogen, 2,000hours lifetime
Power source	AC220V~240V 50/60Hz, AC110V
Power consumption	50VA
Auto -OFF funtion	5~10minutes after
Dimensions	Main body: 265(H) x 230(W) x 360(D)mm
	Remote controller: 20(H) x 64(W)x195(D)mm
Net weight	Main body: 6.0kg
	Remote controller: 180g
Interface	8-pin connector for simultaneous control of chart,mask and R/G filter by a motorized refractor's control unit.
Accessories	1 each of dust cover 1, Polarized metal screen(400mm X 350mm)
	1 polarized glasses, 2 spare fuses and 2 batteries, 1 spare lamp.
Optional accessories	Wall bracket, red/green glasses

• Specification and design are subject to change without prior notice for improvement



unicos
(주)유니코스
Address : 210-4 Maero-Dong, Seo-Gu, Daejeon, Korea
Tel : +82-42-581-0047, 0048
Fax : +82-42-581-0053
http://www.e-unicos.com
E-mail: master@e-unicos.com

Distributed by

